bluephase®

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English

Factory settings

If the enclosed light probe is used, the light intensity is calibrated to the following values:

	bluephase® C8	bluephase® style	bluephase®	bluephase® 20i
Calibrated light intensity	800 mW/cm ² ± 10%	1.100 mW/cm ² ± 10%	1.200 mW/cm ² ± 10%	2.000 – 2.200 mW/cm ²
10 mm light probe, black 636240	—	■	—	—
10 mm light probe, black 608537	■	—	■	—
10>8 mm light probe, black 627389	—	—	—	■
6>2 mm (Pin-Point), black 636241	—	✓	—	—
6>2 mm (Pin-Point), black 608538	✓	—	✓	✓

■ Enclosed in delivery form ✓ Available as accessory

The use of a light probe other than the one provided has a direct influence on the light intensity emitted.

In **parallel-walled** light probes (10 mm), the diameter is equal at both ends. In **focussing** light probes (10>8 mm light probe, Pin-Point light probe 6>2 mm), the diameter of the rear end is larger than that of the light emission window. The incident blue light is thus bundled to a smaller surface, which increases the light intensity emitted.

Pin-Point light probes are suitable for the polymerization of confined areas, such as the attachment of veneers prior to excess removal. For thorough curing, it is necessary to change the light probe.


Important facts: *The larger the distance between the light emission window and the composite restoration, the larger is the surface that is irradiated. The energy emitted by the curing light is spread over a larger surface, but the light intensity decreases accordingly.*

Particularly in the case of posterior restorations, the distance between the cusp and the bottom of the cavity may already be 6 mm or more. Therefore, it is all the more important to always polymerize with the light probe placed directly on the tooth surface or for a longer time.

Double the curing time as follows at the following distances:

Parallel-walled 10 mm standard light probe:
at approx. 9 mm at the latest

10>8 mm light probe:
at approx. 8 mm at the latest

Routine situation: When a composite layer at the bottom of a cavity is cured, the light probe cannot be placed directly on the material.

a : Indicated light intensity at the light emission window (e.g. 1,100 mW/cm²)

b : Decreasing light intensity with increasing distance (e.g. 800 mW/cm²)

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